



ORDER FORM

PYRO LIGHT OPTIMA XL

Light weight wheelchair

Company:

Basic price: **Pyro Light Optima XL without drum brake** € 1.399,-

Basic price: **Pyro Light Optima XL with drum brake** € 1.598,-

Customer:

Max. capacity 170 kg

Date:

▼ Please mark!

Seat	€
Seat width in mm 520 550 580	<input type="radio"/>
Seat height front in mm 500	<input type="radio"/>
400 included front wheel 6", rear wheel 22"	<input type="radio"/> 45,50
425 450 475	<input type="radio"/> 45,50
Seat depth 450 mm (Standard) adjustable from 450 mm up to 500 mm in 25 mm steps	<input checked="" type="radio"/>
Seat cushion Seat width 520 mm	34,-
Seat cushion Seat width 550 mm	40,-
Seat cushion Seat width 580 mm	44,-
Seat upholstery, Nylon	<input checked="" type="radio"/>
Back	
Back support height 420 mm	<input checked="" type="radio"/>
Back support height adjustable 400 - 460 mm	<input checked="" type="radio"/>
Back support with lumbar curve	<input checked="" type="radio"/>
Back support adjustable by strap-system	<input checked="" type="radio"/>
Head support incl. fixation	188,-
Handgrip elevating, removable	<input checked="" type="radio"/>
Back support upholstery, Nylon	<input checked="" type="radio"/>
Sidepanels	
Sidepanels adjustable in depth and height	<input checked="" type="radio"/>
Leg support	
Swing-away leg support with individual footplates plastic, removable, adjustable in height	<input checked="" type="radio"/>
Leg support elongation 100 mm Pair	60,80
Lower leg support elevating left right	107,60
Amputee support left right	157,-
Footplates adjustable in angle and in width	<input checked="" type="radio"/>

* for sw 550 mm additional telescopic tubes are required € 196,-

▼ Please mark!

Tyres	€
Rear wheels 24" x 1 3/8" PU incl. quick-release axle	<input checked="" type="radio"/>
Front wheel 8" 200 x 50 mm, PU	<input checked="" type="radio"/>
Hand rim rubber-coated 24" Pair	110,-
Spoke protection, transparent 24"	145,-
Wheel base extension	72,-
One arm drive only without drum brake left right	325,-
Wheel locks	
Knee lever wheel lock push to lock	<input checked="" type="radio"/>
Wheel lock lever extension	23,-
Knee lever wheel lock for one-arm drive left right	88,60
Accessories	
Togo power push aid incl. mounting set XL version sw max. 550 mm*	2.446,60
Mounting aid for adaptation	189,50
Tilting aid	<input checked="" type="radio"/>
Anti-tip device Pair	87,50
Safety belt XL	30,90
Cane holder	23,70
Therapy tray	215,-
Oxygen/infusion holder and bracket	169,-
APINO Citybag	60,42
APINO Shopper	18,20
Passive light	<input checked="" type="radio"/>
Frame colour	
Silver	<input checked="" type="radio"/>

The specified dimensions are subject to manufacturing tolerances of +/- 10 mm

Basic equipment

Optional, no extra charge